

IMPORTANT! If your bed is not operating properly, you may need to reset it with the remote.

Resetting Your SleepSafe® Bed



Resetting a TWIN SleepSafe® Bed

1. Locate the DOWN arrow button for the Hi-Lo function of the Remote Control.
2. Press button until you hear the bed click once, and continue to press button for 6 to 8 seconds until you hear the bed click a second time.



Resetting a FULL or QUEEN SleepSafe® Bed

1. Locate the UP and DOWN Arrow button for the Hi-Lo function of the Remote Control.
2. Press both buttons at the same time until you hear the bed click once.

These procedures should reset the bed's T-MOTION® control box. If your bed is not operating, please call us for technical support.



SleepSafeBed.com / 866-852-2337

Made in USA