

IMPORTANT! When ordering your bed, please indicate with or without padding. If padded, see options below.

SleepSafe® Bed Padding Options



SleepSafe II® - MEDIUM BED
Light Blue padding over mesh screened windows.



SleepSafe II® - MEDIUM BED
Light Blue padding around mesh screened windows.

Other Color Choices



Gray



Cream



Taupe



Gray



Cream



Taupe

Materials:

- 70 Denier Nylon Fabric Cover with a Fluid Resistant Healthcare Coating
- 3/4" High Density Foam Padding

Padding features:

- Bacteria resistant
- Fungal resistant
- Non-allergenic
- Malodor Protection
- Flame Retardant
- Stain Resistant
- Fluid Resistant
- Tear Resistant
- Anti-Static

The safety rail pads extend up and over the rails and attach to the outside of the safety rails. The headboard and footboard pads extend to the top and to the edges of the legs of headboard and footboard.

Padding questions?


- Please contact customer service for any specific information about padding.

Padding attachment:

- Industrial strength Velcro is used to secure the padding to inside panels and safety rails.

Padding maintenance:

- For cleaning, use warm water and a mild cleanser. If the padding is over the windows, it can be removed and cleaned. If the padding is around the windows, please do not remove. It is best to clean this padding in its place.



Important notes:

Padding is installed at the factory and will cover ALL 4 sides. Padding only 1 or 2 sides is NOT an option. Padding removal will cause a gap between the mattress and the bed's interior.



SleepSafeBed.com / 866-852-2337

Made in USA / Call for DME information / FDA Registered