

IMPORTANT! When ordering your bed, please indicate with or without padding. If padded, see options below.

SleepSafe® Bed Padding Options



Padding Color Choices



Materials:

- 70 Denier Nylon Fabric Cover with a Fluid Resistant Healthcare Coating
- 3/4" High Density Foam Padding

Padding features:

- Bacteria resistant
- Fungal resistant
- Non-allergenic
- Malodor Protection
- Flame Retardant
- Stain Resistant
- Fluid Resistant
- Tear Resistant
- Anti-Static
- Padding extends to the top of the headboard and footboard.
- Padding extends up and over the safety rails and attaches to the outside.
- Padding covers the corner posts.

Padding questions?


- Please contact customer service for any specific information about padding.

Padding attachment:

- Industrial strength Velcro® is used to secure the padding to inside panels and safety rails.

Padding maintenance:

- For cleaning, use warm water and a mild cleanser. If the padding is over the windows, it can be removed and cleaned. If the padding is around the windows, please do not remove. It is best to clean this padding in its place.



Important notes:

Padding is installed at the factory and will cover ALL 4 sides. Padding only 1 or 2 sides is NOT an option. Padding removal will cause a gap between the mattress and the bed's interior.



SleepSafeBed.com / 866-852-2337

Made in USA / Call for DME information / FDA Registered