



SleepSafe Beds, LLC
3629 Reed Creek Drive
Bassett, VA 24055

Toll Free: 866.852.2337
Local: 276.627.0088
Fax: 276.627.0234

SleepSafeBed.com

IMPORTANT MESSAGE WHEN ORDERING PADDING FOR A SLEEPSAFE® BED

Please read this section carefully

- **IMPORTANT:** Please specify padding with your initial *SleepSafe® Bed* order.
- When ordering padding, please indicate if you want the padding to cover the windows in the safety rails, or if you want the padding to be installed around the windows in the safety rails. Alternatively, the option of solid safety rails is available.
- A *SleepSafe® Bed* ordered with pads is fitted with a mattress to accommodate the padding.
- A *SleepSafe® Bed* ordered without pads is fitted with a wider mattress.
- If mesh windows are requested, in lieu of plexi-glass, padding will automatically be installed around the windows.
- If you want to retrofit padding to an existing *SleepSafe® Bed*, a conversion kit is available. Contact customer service at 866.852.2337 for details.



PADDING INSTALLED
OVER WINDOWS



PADDING INSTALLED
AROUND WINDOWS



Other Notes:

- If you order padding for your *SleepSafe® Bed*, it will be installed at the factory and will cover ALL FOUR sides. Padding only one or two sides is NOT an option.
- If you order padding for your *SleepSafe® Bed*, it is very important to NEVER use the bed if you have removed one or more of the pads. If you remove ANY of the pads, a gap will be created between the mattress and the sides of the bed potentially creating a risk for entrapment.